

# Stellar M22

Shine on




## ResurFX™

For True Fractional Non-Ablative Skin Resurfacing


Needs just one pass for effective skin resurfacing, saving time and protecting patient's skin

- ◆ Extremely versatile state-of-the-art CoolScan™ scanner, for customizable treatments
- ◆ Unique algorithm for controlled non-sequential pulsing ensures precise, homogenous and safe energy delivery for reproducible results
- ◆ Efficacious and fast, with up to 20% coverage in one pass
- ◆ Cost effective upgrade with NO disposables



# Fractional non-ablative resurfacing treatments in a single pass

## Typical Clinical Results



Courtesy of Matteo Tretti Clementoni, MD.



Courtesy of Matteo Tretti Clementoni, MD.



Courtesy of Mitchel P. Goldman, MD



Courtesy of Arielle N.B. Kauvar, MD



Courtesy of Thierry Fusade, MD



Courtesy of Olga Yarish, MD

**Warnings and risks:** Stellar M22™ ResurfXTM is contraindicated for patients with active infection, viral, fungal or bacterial diseases. The use of ResurfXTM could cause redness, swelling, scarring, change of pigmentation. See the system user manual for a complete list of contraindications and risks.

### ResurfXTM Specifications

Wavelength	1565 nm
Energy	Up to 70 mJ / microbeam
Scan Size	Up to 18 mm diameter
Tips	Sapphire tip; Precision tip

Scan Shapes	Horizontal line, vertical line, square, horizontal rectangle, vertical rectangle, circle, donut, hexagon
Beam Density	Up to 500 microbeams / cm <sup>2</sup>
Repetition Rate	Up to 2 Hz
Cooling	Continuous contact cooling

© 2020 Lumenis Ltd., All Rights Reserved. The Lumenis Group of Companies, PB-00002780 Rev A. Lumenis, its logo, Stellar M22, ResurfFX and CoolScan are trademarks or registered trademarks of the Lumenis Group of Companies. Specifications are subject to change without notice Lumenis is ISO 13485:2016 certified.

**Lumenis Ltd.**  
Yokneam Industrial Park  
6 Hakidma Street P.O.B. 240  
Yokneam 2069204, Israel  
Tel: +972-4-959-9000

**EMEA**  
Dreieich Dreieichenhain,  
Germany  
T +49 6103 8335 0  
F +49 6103 8335 300

**JAPAN**  
Tokyo, Japan  
T +81 3 6743 8300  
F +81 3 6743 8301

**ASIA PACIFIC**  
Gurgaon, India  
T +91 124 422 07 95  
Kowloon, Hong Kong  
T +852 217 428 00  
F +852 272 251 51

**Lumenis (Germany) GmbH**  
Heinrich-Hertz-Str 3 D-63303  
Dreieich-Dreieichenhain  
Germany  
Tel: +49 (0) 6103 8335 0

**Roma, Italy**  
T +39 06 90 75 230  
F +39 06 90 75 269

**AMERICAS**  
San Jose, CA, USA  
T +1 408 764 3000  
+1 877 586 3647  
F +1 408 764 3999

**Hertfordshire, UK**  
T +44 20 8736 4110  
F +44 20 8736 4119

**Sydney, Australia**  
T +61 1800 5863 64



Rx ONLY

© 2020 LUMENIS LTD., ALL RIGHTS RESERVED



LUMENIS.COM/Stellar-M22

PB-00002780 Rev A